

浙江中在医疗科技有限公司

通用移液吸头兼容移液枪品牌适配表

品牌	规格	量程	P0010	P0200	P0300	P1000	P1250	PX0010	PX0100	PX0200	PX0500	PX1000	PX1250	备注	
美国 赛默飞/thermo	Single Channel	F3系列 1-10u1	√					√							
		F3系列 2-20u1		√					√	√					
		F3系列 5-50u1		√						√					
		F3系列 10-100u1		√						√					
		F3系列 20-200u1		√						√	√				
		F3系列 100-1000u1					√	√				√	√	√	
		F2系列 0.5-5u1	√												
	F2系列 1-10u1	√						√							
Multi Channel	F2系列 排枪1-10u1	√						√							
赛多利斯/sartorius	Single Channel	proline系列0.1-2.5u1													
		proline系列0.5-10u1	√						√						
		proline系列2-20u1		√						√					
		proline系列10-100u1		√						√					
		proline系列20-200u1		√	√					√	√				
	proline系列100-1000u1					√	√					√	√	需稍用力 插入吸头	
Multi Channel	排枪0.5-10u1	√						√							
普兰德/brand	Single Channel	0.1-2.5u1	√					√							
		1-10u1	√					√							
		2-20u1		√						√					
		10-100u1		√											
		20-200u1		√	√					√	√				
		100-1000u1					√	√				√	√	√	

品牌	规格	量程	P0010	P0200	P0300	P1000	P1250	PX0010	PX0100	PX0200	PX0500	PX1000	PX1250	备注	
吉尔森/gilson	Single Channel	PN型2-20u1		√					√						
		PN型10-100u1		√					√						
		PN型20-200u1			√						√				
		PN型100-1000u1				√	√					√	√	√	
艾本德/eppendorf	Single Channel	0.1-2.5u1	√					√							
		1-10u1	√					√							
		2-20u1		√						√					
		10-100u1		√						√					
		20-200u1		√	√						√				
	100-1000u1				√	√					√	√	√		
	Multi Channel	0.5-10u1	√						√						
		5-100u1		√						√					
		50-100u1		√						√					
		20-300u1		√	√						√				
30-300u1			√	√											
50-1200				√	√										
大龙/dragonlab	Single Channel	0.1-2.5u1	√					√							
		1-10u1	√					√							
		2-20u1		√						√					
		10-100u1		√						√					
		20-200u1		√	√						√				
		100-1000u1				√	√					√	√	√	